


BUS/SCS nurse call unit

signalling units for nurses' station

Antimicrobial⁽¹⁾



Pack	Cat.Nos	Antimicrobial ⁽¹⁾ control units
		<p>Call and presence display for all bedrooms and rooms in order of priority (max. 150)</p> <p>4-module wiring mechanism supplied with support</p> <p>Installation in 50 mm deep flush-mounting box or on inclined plane Cat.No 0 766 12/13 with alphanumeric display and adjustable audible signal</p> <p>Scrolled display of the calls and nurse visits in progress</p> <p>Fitted with a malfunction indicator</p> <p>Used in conjunction with the interphone unit</p> <p>Cat.No 0 766 08 allows:</p> <ul style="list-style-type: none"> - communication with the bedrooms issuing calls and with a nurse present - call acknowledgement at communication end - calls to be made to all rooms with a nurse present <p>Connects directly to the BUS/SCS cable</p> <p>Cat. Nos 0 492 32/33</p> <p>Main control unit</p> <p>1 0 766 11 For installation in the nurses' station</p> <p>Secondary control unit</p> <p>1 0 766 09 Can be used in addition to the main control unit in another room (break room, dining room, etc.) or in a corridor</p>
		<p>BUS/SCS cables</p> <p>1 0 492 33 Length 200 m (reel)</p> <p>Non-halogen cable</p>
		<p>Fixing accessories</p> <p>For control units</p> <p>1 0 766 12 Inclined plane for fixing control unit only onto table</p> <p>1 0 766 13 Inclined plane for fixing control unit + interphone unit onto table</p> <p>For control units, door units and interphone units</p> <p>1 0 766 14 Surface mounting frame for refurbishment</p>
		<p>Call transfer interface to DECT</p> <p>1 0 766 19 Interface for forwarding information as ESPA 444 protocol to the DECT system (text messages only)</p> <p>2 x 17.5 mm DIN modules</p>

Pack	Cat.Nos	Traceability on PC
1	0 766 18	<p>Traceability software</p> <p>Memorises and categorises events by information type (calls, nurse presence, acknowledgement, etc.)</p> <p>Time stamping of events</p> <p>Export of the elements as .csv files which can be used on a PC</p> <p>Works with traceability interface Cat.No 0 766 17 to be ordered separately</p>
1	0 766 17	<p>Traceability interface</p> <p>Works with call traceability software Cat.No 0 766 18</p> <p>2 x 17.5 mm DIN modules</p> <p>1: Contains a silver compound to prevent the growth of bacteria on the surface of the plate</p>

Wiring diagram example, see e-catalogue



More information in our "Solutions for retirement homes & healthcare establishments" brochure available on www.docexport.legrand.com

